

ACCESSWALK® Specifications & Warranty

- Dimensions of ACCESSWALK® shipping in a cardboard box: 39" × 11" × 38" for 40 lb. (98 cm x 28 cm x 95.5 cm for 18 kg)
- ACCESSWALK® fits through standard ADA doorways
- Lightweight anodized aluminum frame
- Foldable for storage and transportation
- Folded Up Dimensions: (without 12" front wheels & back bar): 33.5" L (85.1 cm) x 34.45" W (87.5 cm) x 10" H (25.4 cm)
- Folded Up Dimensions: (with 12" front wheels): 33.5" L (85.1 cm) x 34.45" W (87.5 cm) x 17" H (43.2 cm)
- Width of ACCESSWALK®: 34.45" (87.5 cm)
- Length: 30" (76 cm)
- Net Weight: 37.5 lb. (17 kg)
- Gross Weight: 39.7 lb. (18 kg)
- Maximum weight capacity: 285 lb. (130 kg)
- Folding seat and textile storage basket
- Seat Material: Polyester mesh incl. foamed cushion
- Seat Height: 23" (59 cm)
- Seat Dimensions: 16.3" wide x 11.8" deep (41.5 cm x 30 cm)
- Parking brake: Yes, activated with the break hand brake
- Turning Radius Circle: 35.43" (90 cm)
- Front Wheel Stainless Steel casters/axel
- Yellow "SUNNY" Pneumatic rubber wheels: Two front wheels of 12" diameter (30.5 cm) and two back wheels of 10" diameter (25.4 cm)

ACCESSREC

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: sales@AccessRec.com

Web: www.ACCESSRECC.com

- Adjustable/removable handles from 32.5" to 37" high (82.55 cm to 94 cm)
- Width Between Handles: 18" (45.7 cm)
- Handles material: Plastic
- Carry Basket: 11 lb. payload (5 kg)
- Carry Basket Material: Mesh fabric
- Warranty: Limited Warranty of 3 Years

ACCESSREC

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: sales@AccessRec.com

Web: www.ACCESSREcc.com

ACCESSREC Rollator Walker – Limited Warranty

ACCESSREC provides three (3) years limited warranty from the date of purchase, covering defects in materials and workmanship for the rollator frame only.

Warranty Exclusions

The following components are not covered under this warranty:

Seat cover material

Carry basket

Brakes

Plastic handles or other plastic spare parts

Wheels

All warranty claims and limitations will be assessed solely by ACCESSREC.

Warranty Limitations

The three (3) years limited warranty does not cover damage to the rollator frame caused by:

Normal wear and tear

Improper use, maintenance, or storage

Accidents or natural disasters

Unauthorized modifications made without prior written consent from ACCESSREC

Unauthorized repairs performed without prior written approval from ACCESSREC

For any warranty inquiries or claims, please contact ACCESSREC directly.

ACCESSREC

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: sales@AccessRec.com

Web: www.ACCESSRECC.com